

<b>Specimen Type:</b>	<b>Washing</b> <i>(Bladder, Bronchial, Gastrointestinal tract, Pelvic, Peritoneal)</i>
<i>Patient Preparation:</i>	As ordered by M.D.
<i>Specimen Collection:</i>	The collected material should be transferred to a cytology specimen container.
<i>Specimen Requirements:</i>	Submit in saline. Cytology fixative may be added.
<i>Handling/Transport:</i>	If submitted fresh, submit to laboratory immediately, refrigerate if delayed.

---